



Remote Access

There are multiple ways you can access PowerSoftMD from a remote location.

- 1) Leave a computer on at the office and install a remote access program.
Make sure the Windows screen saver is turned on, requiring a password for access.
We also highly recommend the computer be powered using through a battery backup.
There are many services available; currently, the most cost effective one we have found is www.remotepc.com (about \$55/year). RemotePC supports both **Windows** and **Apple** type remote computers as well as iOS and Android devices. Remote PC directly handles remote printing and claims to be HIPAA compliant.

Other services include: LogMeIn, GoToMyPC, and TeamViewer.
- 2) If your main computer is a Windows Server such as Windows 2016 Server or Windows Server 2019, you can have your hardware/network technician set up Remote Desktop Services allowing you to use using **RDP** on your remote computer for access. For remote printing support, we recommend using a program TSPrint from www.terminalworks.com. For remote scanning, we recommend using a program called Scan Redirector from www.incentivespro.com. RDP is available for both **Windows** and **Apple** type computers. Again, you will need your hardware/network technician set up this type of access.
- 3) If you start using a Tablet-PC or iPad you may like to try the Tablet style PowerSoftMD Menu. From the top of the Primary Menu click "**Tablet-Menu**", if you don't like it you can switch back by clicking "Desk Top Menu".
- 4) Using a **VPN** connection [Does Not Work](#). VPN is simply **TOO SLOW**, regardless of your internet connection speeds.